ENGLISH

Technical Specifications

- Power supply: 2 x 1.5 V AAA batteries
- Pulse frequency: 0-1000 Hz
- Pulse width: 20-400 µs
- Maximum power per pulse: ≤ 300 mJ
- Magnetic field strength: 0.01 T

Instructions

- Insert the batteries into the compartment and connect the Electromagnetic Neck & Back Massager Innovagoods® with the front ON/OFF button.
- Press the ON button on the side to switch on the device in automatic massage mode 1. Press the MODE button to select the following massage programs that will appear on the screen with the following display indications:

1. Automatic (1)
2. Acupuncture (2)
3. Massage (3)
4. Shiatsu (1+3)
5. Multiple (1+2)
6. Impacts (1+2+3)

- Intensity adjustable up to 16 levels with the HIGH/LOW buttons on the side.
- Place the Electromagnetic Neck & Back Massager Innovagoods® on your neck and make sure the stainless steel tipping pads are in contact with your skin.
- Put the silicone patches on the electrodes, connect the cable to the SOCKET plug and stick the patches on the area you want to massage: back, shoulders, lower back, etc.
- You can use the Electromagnetic Neck & Back Massager Innovagoods on the neck or simultaneously on the neck and back by installing the electrodes.
- Recommended time of use: 10-20 minutes daily.

WARNINGS

If you feel any discomfort when using the device, discontinue use and consult a doctor. Do not use the device if you have heart problems, epilepsy or if you are pregnant. Do not use the device on areas where the skin is infected or sore. Do not mix new and used batteries in the compartment. If you are not going to use the unit for a long period of time, remove the battery. Do not expose batteries to high temperature or direct sunlight. Do not attempt to open them. do not dispose of them in fire, there is risk of explosion. Clean the device with a slightly damp cloth. Do not use abrasive products or materials to clean it. Do not immerse in liquid. Do not handle it with wet hands. Do not use the device close to a source of water or in damp places. Suitable for domestic use only. Use the device for its intended purpose only. Store the product in a dry place with no dust. To avoid risks, do not use it if any parts are damaged and do not attempt to disassemble or repair it. This device is not designed for use by children or people with physical and mental disabilities, unless they are fully supervised by an adult who is responsible for their safety. This product is not household waste, so follow the environmental standards for recycling. Discard the product in accordance with the European Directive 2002/96/EC. Enquire about your local recycling legislation for electrical and electronic products and follow these steps. Do not dispose of this product with household waste. Following these instructions will help prevent harmful consequences to personal and environmental safety.

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